



BETAtex
Medical mask
disposable, with earloops
TD-30

Attachment No. 3
point 2 pp. 2.3

Revision I
2020.07.31

Page 1/1

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- Protects operative field against exfoliated epidermis, hair and droplet infections

FEATURES:

- Made of three layers of non-woven
- Modelling nose-piece
- In accordance with the standard EN 14683 type IIR
- Splash resistant
- BFE \geq 98%
- Dimensions: 175 mm x 95 mm
- With earloops
- Packaging: ergonomic cardboard box to remove individual pieces
- Colours: green, blue, white, pink
- Latex-free
- Single-use
- Non-sterile



SIZE AND REFERENCE

| CATALOGUE NUMBER | COLOUR | INDIVIDUAL PACKAGING | INTERMEDIATE PACKAGING | BULK PACKAGING |
|------------------|--------|----------------------|------------------------|----------------|
| BT-046-R31G | green | box | 50 pcs | 20 x 50 pcs |
| BT-046-R31B | blue | box | 50 pcs | 20 x 50 pcs |
| BT-046-R31W | white | box | 50 pcs | 20 x 50 pcs |
| BT-046-R31P | pink | box | 50 pcs | 20 x 50 pcs |